

BACTERIAL REVERSE MUTATION TEST (AMES)

The **Bacterial Reverse Mutation** test assesses the mutagenic potential of a compound and is commonly employed as an initial **screen for genotoxic activity** and, in particular, for mutation-inducing activity, which involve substitution, addition or deletion of one of the DNA base pairs.

STRAINS

Ames testing uses several strains of the bacterium *Salmonella typhimurium* which carry a defective (mutant) gene that renders them unable to synthesize the amino acid histidine. The Ames test investigates the potential of the test compound to result in a **back mutation** that causes the gene to regain its function and grow in a histidine-free medium.

At least five strains of bacteria should be used :

- **Base-pair substitution**: *S. typhimurium* TA1535 and TA100 and TA102 or *E. coli* WP2 uvrA or *E. coli* WP2 uvrA (pKM101)
- **Frame-shift**: *S. typhimurium* TA1537 or TA97 or TA97a and TA98: frame-shift

METABOLIZATION

Mutagenic potential can be investigated in the Ames test in the presence or absence of a **metabolizing system** (rat liver S9 fraction) to identify pro-mutagens as well as directly acting mutagens.

CONDITIONS

The standard methods apply for most pure chemicals and consist of **incorporation method, preincubation method, fluctuation method and suspension method**. An additional method referred as **Treat and Wash method**, applies for proteins natural or synthesized, mixture, plant extract, antibiotics, or cytotoxic compound (oncologic).



Optimized genotoxicity characterization for limited sample quantities

For **non GLP exploratory screening** *S. Typhimurium* TA98 (frameshift mutation) and TA100 (base-pair substitution) are two common strains used for this reduced Ames test since they both have *rfa* mutations (a defective lipopolysaccharide layer that makes bacteria more permeable to larger molecules), *uvrB* mutations (elimination of excision repair of DNA damage) and pKM101 plasmid (increasing error-prone repair of DNA damage).

GenEvolutionN has developed **High Throughput miniaturized conditions (non GLP)** when a limited amount of sample is available (discovery genotoxicity screening and impurity genotoxicity characterization), such as **MINI Ames test or NANO® Ames test**.

| | |
|---------------------------|--|
| Test model | <i>Salmonella typhimurium</i> TA 98, TA100, TA1535, TA1537, TA102 (optional <i>E.coli</i> WP2 UvrA) |
| Assay Controls | Negative control: Aqueous and organic solvents Positive control: Strain-specific positive control |
| Exposure time | 48 – 72 hours |
| Test item quantity | Standard condition (incorporation or preincubation): 1 g Treat & Wash conditions: 250 mg MINI Ames, 25 wells plate (agar): 250 mg MINI Ames, MPF 384 wells plate (liquid): 65 mg NANO® Ames (25 and 96 wells plate): 35 µg , also applicable for unknown impurities in solution (no additional synthesis of compound) |
| Endpoint | Growth of revertant colonies combined with genomic sequencing for NANO® Ames |
| Data delivery | Mean number of revertant colonies per plate <ul style="list-style-type: none">• 2-fold increase with vehicle control (TA 98, TA 100, TA97, TA97a, TA102, and WP2 UvrA)• 3-fold increase with vehicle control (TA 1535 and TA 1537) |
| Timeline | 2 – 3 weeks |
| Regulatory statut | GLP compliance for Standard and Treat & Wash conditions with Formulation Analysis (optional) |